





YOU CAN DO IT

WITHOUT LANCETS¹

Moving from the
FreeStyle Libre system to
the FreeStyle Libre 2 system?



We're pleased to introduce you to the second generation of our sensor-based technology for people living with diabetes. The FreeStyle Libre 2 system comes with the same great features of the FreeStyle Libre system and more.


FreeStyle
Libre 2
FLASH GLUCOSE MONITORING SYSTEM


life. to the fullest.®
Abbott

Images are for illustrative purposes only. Not actual patient or data.
1. Scanning the sensor does not require lancets.

Moving to the FreeStyle Libre 2 system is easy.
Follow the guidance below.

**TALK TO YOUR DOCTOR AND CONFIRM YOUR PRESCRIPTION
HAS BEEN CHANGED TO FREESTYLE LIBRE 2 SENSORS**



**CHOOSE WHICH DEVICE YOU WOULD LIKE TO RECEIVE
ALARMS ON?¹**



PHONE



YOU'RE READY TO GO!

Using the FreeStyle LibreLink app, you can start scanning FreeStyle Libre 2 sensors straight away²

Please check that their FreeStyle LibreLink app version is 2.5 or above. See compatibility guide for details.



READER



REQUEST A READER

Go to the Reader Replacement Portal and complete the online request form

[FreeStyleLibre.co.uk/replacement](https://www.freestylere.com/FreeStyleLibre.co.uk/replacement)



REQUEST VALIDATION

- Enter your current FreeStyle Libre reader serial number
- No requirement to send the current reader back



**DELIVERED DIRECTLY
TO YOU**

A replacement reader is delivered directly to you as the FreeStyle Libre user



Images are for illustrative purposes only. Not actual patient data.

1. Patients choose which device they want to receive alarms: FreeStyle Libre 2 reader or FreeStyle LibreLink app. They must start their FreeStyle Libre 2 sensor with that selected device. Once the patient scans their FreeStyle Libre 2 sensor with that device, they can receive alarms only on that device. **2.** The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink requires registration with LibreView.

EDUCATION AND SUPPORT FROM ABBOTT – SUPPORTING USERS THROUGH THE MOVE



FreeStyle Academy



Webinars



Website information



New starter guide

Prescribing information.



THE FREESTYLE LIBRE 2 SENSOR



Sensor (1 sensor lasts up to 14 days)
PIP Code: 416-3416

Recommended quantity

Discover more support and guidance for you.
Visit us at: **www.FreeStyleDiabetes.co.uk**

What is the difference between the FreeStyle Libre system and the FreeStyle Libre 2 system?

Just like the FreeStyle Libre system, the FreeStyle Libre 2 system is discreet, convenient and easy to use¹, now with added benefits:



No finger pricks²

There's no need to prick your finger when using the FreeStyle Libre 2 system - even when glucose is low, falling or rapidly changing.



Excellent accuracy^{1,3}

Outstanding accuracy in the low glucose range when it matters most



Optional alarms

Peace of mind when it's needed most, with three optimal real-time alarms for those who need them.

- Low glucose alarm
- High glucose alarm
- Signal loss alarm


FreeStyle
Libre 2
FLASH GLUCOSE MONITORING SYSTEM


life. to the fullest.®
Abbott

Images are for illustrative purposes only.

1. Data on file, Abbott Diabetes Care, Inc. **2.** Finger pricks are required if your glucose readings and alarms do not match symptoms or expectations. **3.** Alva S, et al. Accuracy of a 14-Day Factory-Calibrated Continuous Glucose Monitoring System With Advanced Algorithm in Pediatric and Adult Population With Diabetes. Journal of Diabetes Science and Technology. September 2020. doi:10.1177/1932296820958754.

© 2020 Abbott. FreeStyle, Libre, and related brand marks are marks of Abbott. ADC-32039 v1.0 12/20.